

Dr's Light Clever

Incomparable Dual Wave Light Quality

Dr's Light Clever is a new conceptual curing light applied many advantages for user's demands.
Just press the button and experience new light technology with Dr's Light Clever.




Incomparable Dual Wave Light Quality

Dr's Light Clever is a new conceptual curing light applied many advantages for user's demands.
Just press the button and experience new light technology with Dr's Light Clever.



330° head twist and slim head



Autoclaveable and interchangeable Guide-cap




Interchangeable battery pack

1 Features of Dr's Light Clever Dual


- Dual wavelength LED**
Usable all kinds of materials in the field with 405nm & 460nm peak wavelength LED.
- Focused light**
Increase efficiency of light.
- Autoclaveable and interchangeable Guide-cap**
Guide-cap is autoclaveable and interchangeable to prevent cross infection.
- 4 curing light modes**
User can select High, Low, Softstart, Pulse mode as user wants.
- Time setting available (Softstart is fixed for 20 sec.)**
User can select time setting as needs.
- 330° head twist and slim head**
It makes easy to get approach into the oral.
- Angled handpiece design**
Most familiar design for user.
- OLED screen**
It shows device states on the OLED screen with clearly.
- Interchangeable battery pack**
User can interchange the battery easy.
(Doesn't need to send device to factory to change the battery)
- Big capacity battery**
Big capacity battery makes long time use.

2 Comparison with Other products



3 Various curing modes


Selectable light intensity and time as you want




4 Specification

Model : CL-DC20 / Wave length : 400~490(nm) / Guide-cap : Autoclaveable /
Power source : Lithium-ion-battery / Light source : High Power Dual wavelength LED /
Weight (Handpiece) : Approx. 130g

5 Accessory



Guide-cap (x 3 pcs)



Battery