

# C-SAILOR PRO Dental Implantation System

Intuitive  
Powerful  
Safe & Reliable

03 Dental Implantation Systems



## Patent design, delicate contracted

This product has been applied layout-design patent, Counterfeit must be investigated.



The large size screen ensures that the treatment parameters are visible at any time

Model: C-Sailor

Dental Implantation Systems

Order code: C-SAILOR PRO

### Technical:

Maximum power:	130VA
Power supply:	AC110V 60Hz / 230V 50Hz
Motor speed:	300--40000rpm
Adjustable torque:	5Ncm--80Ncm (20:1)
Water supply volume:	4 level adjustable, up to 150ml/min
Light source illumination:	30000lux
Motor:	surgical brushless motor
Volume(cm):	42x38x26cm
Weight(kg)/PCS	5.5kg
Packing material	Aluminum box
QTY per master carton	2 pc/ctn



Fiber optic contra angle  
High brightness LED light  
source up to 30000lux




Built-in water supply system  
4 shifts water supply  
up to 150ml/min

04 Dental Implantation Systems



-Full touch 7 inch screen-  
Large size screen ensures that  
the treatment parameters are  
visible at any time



Full function pedal controller  
Aseptic operation can be achieved  
giving the maximum freedom of the handle



# C-SAILOR Dental Implantation System

Efficient  
Big screen  
Stable & Safe



The first Chinese brand implant machine which get license from Chinese government for producing and selling

05 Dental Implantation Systems



This product has been applied layout-design patent,  
Counterfeit must be investigated:NO.ZL201230103509.0

Model: C-Sailor

## Dental Implantation Systems

Order code: C-SAILOR

### Technical:

Maximum output power	120VA
Power	AC110V 60Hz / 230V 50Hz
Adjustable torque,range	5.0-55N.cm(20:1)
BLDC Speed	300rpm-5000rpm
Infusion flow(maximum)	150ml/min
Motor	Surgical brushless motor
Volume(cm)	40X33X26.5
Weight(kg)/PCS	8.2
Packing material	Aluminum box
QTY per master carton	2 pc/ctn



High performance control system  
10 kinds of the reduction contra angle to choose  
Built-in 10 setting program



06 Dental Implantation Systems

### Big Screen

Clear and easy to understand



- Coolant control
- Program control
- Forward/Reverse



### Multi functional foot pedal

Could achieve all the operation function  
Easy and flexible



### Built-in water supply system

5 shifts water supply  
up to 150ml/min



Model: C-TW1

## Dental Implant Torque Wrench



Order code: C-TW1 Kit

### Feature:


1. 10/15/20/25/30/35N torque gears
2. 16 kinds of screwdrivers
3. Used and upgraded the COXO implanting contra angle

### Technical:

Torque: 10-35N.cm

Screwdriver type:  $\phi$ 2.35

Volume(cm) 13.8x13.8x15



Order code: C-TW1

The conventional design of COXO implanting contra angle has been used with better quality and longer service life



Order code: C-SK1

16 kinds of screwdrivers of different specifications fully compatible with more than 30 mainstream implant systems

07  
Dental Implant  
Torque Wrench

Model: C-TW2

## Dental Implant Torque Wrench



Order code: C-TW2 Kit

### Feature:

1. Equipped with 16 different screwdrivers,
2. More optimized structural design, more in line with ergonomics and the actual needs of clinicians.
3. CNC precision machining to ensure accurate torque.
4. Aluminum alloy box, can be sterilized at high temperature and high pressure
5. Paired with silicone pins of different colors

### Technical:

Torque: 15-70N.cm

Screwdriver type:  $\phi$ 7.5

Volume(cm) 15.7x12.5x5



Order code: C-TW2

More optimized structural design, more in line with ergonomics and the actual needs of clinicians



Order code: C-SK2

16 kinds of screwdrivers of different specifications, fully compatible with more than 30 mainstream implant systems

# Endo Series





**Head Office:** Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE

**Phone:** +971 4 569 3474

**Mobile:** +971 50 912 4577

**Email:** [admin@skyline-med.com](mailto:admin@skyline-med.com)

