



SKYLINE


MEDICAL REQUISITES TRADING LLC

360°CRYOLIPOLYSIS CRYOCOOL

The Newest Generation Applicator



5 APPLICATOR FULL BODY SHAPING



Treatment Head P1



Treatment Head P2



5 APPLICATOR

MORE TREATMENT AREAS MORE EFFECTIVE

Treatment Head P3



Treatment Head P4





Treatment Head P5



ADVANTAGES -12° SUPPER COOLING





Cooling on both sides



Traditional Cryolipolysis by two Chilling Plates, only 40% effective for treatment area

360° cooling



NEW 360° cooling applicator
100% effective for treatment area
The lowest temperature can reach -12 degrees

WORKING MECHANISM


After the adipose tissue of the target care site is sucked into the cooling probe, the cooling probe releases a 360-degree surrounding cooling energy to act on the fat cells without damaging any surrounding tissue. Cryocool's steady and powerful cooling energy triggers a natural cell death process called "apoptosis", after which apoptotic fat cells are excreted through the body's natural metabolic process. Cryocool targets fat cells for permanent removal from the body, resulting in significant fat loss in the treated area.




Many of us have stubborn fat despite diet and exercise.



Warm-up technology to accelerate blood circulation



Applicator 360 degrees around fat cells without damaging any surrounding tissue




In the weeks to follow, your body naturally process the fat and eliminates these dead cells.




Cryocool procedure results are long-term, as treated fat cells are gone for good.

BEFORE & AFTER


0day



30days



90days



Gender: Female

Age: 32

Treatment area: abdomen

Treatments times: 3



Gender: Female

Age: 35

Treatment area: abdomen

Treatments times: 2



Gender: Female

Age: 29

Treatment area: thigh

Treatments times: 4

5-1 APPLICATOR



Abdomen/Side

Back waist

Back of thigh



5-2 APPLICATOR



■ Abdomen

■ Back

5-3 APPLICATOR






- Upper Arm
- Abdomen/Side
- Back
- Thigh



5-4 APPLICATOR



- Thigh(front, back, left and right)
- Waist side



5-5 APPLICATOR



- Double Chin
- Bra Fat
- inside of knee



ADVANTAGES




■ **Medical silicone, safe and comfortable**



■ **Medical filter cotton, effectively filter debris**

ADVANTAGES HUMANIZED TOUCH SCREEN




UI INTERFACE

Smart mode



Fast mode



Select different size applicator

Select cooling mode

Choose fat thickness


Set treatment time

Start/Stop

Adjust vacuum level

Adjust temperature
(cooling/heating)

ADVANTAGES FOLDABLE SCREEN



SPECIFICATION


360 CRYOLIPOLYSIS CRYOCOOL

Input power:	1200W
Frozen Head Pressure:	0-60Kpa
Frozen Head Screen:	-12°C -1°C
Treatment Handle's Type:	5 applicator in different size
Machine Screen:	15.6 inch Collapse screen, 30°-73° adjustable.
Cooling System:	Semiconductor cooling + water cooling +air cooling
Dimension:	41*58*154cm




BEFORE & AFTER


0day



30days



90days



Gender: Female
Age: 32
Treatment area: abdomen
Treatments times: 3



Gender: Female
Age: 35
Treatment area: abdomen
Treatments times: 2



Gender: Female
Age: 29
Treatment area: thigh
Treatments times: 4

360° CRYOLIPOLYSIS CRYOCOOL





Head Office: Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE
Phone: +971 4 569 3474 **Mobile:** +971 50 912 4577 **Email:** info@skyline-med.com

