

# C-SMART-I PRO

Endo motor  
With LED light  
With Apex locator



## Big color LCD screen

Give you what you want to know



## Touch button design

Simplicity and easy to use



## Built-in mainstream system

Not only have 6 mainstream brands  
But also have 7 memory program setting.



## High capacity battery

- Can work continuously for a long time
- Rechargeable



150~600rpm

## Impored Special Motor

Mute·Big torque



0.6~5.2N.cm

## Motor with Apex locator

Three working mode


M1: Apex locator mode

M2: Endo motor mode

M3: Dual mode

## Mini contra angle 16:1

Smaller head, more easy to operate



## LED lighting system;

Detachable design, powerful LED,  
give you a clear vision

