



SKYLINE
MEDICAL REQUISITES TRADING LLC

LED HANDPIECE WITH GENERATOR

Model: CX207-F

New upgrade, new LED bulb, new integrated head



Self-luminous

Graceful body,simple,protable
easy clean,easy control

Balanced impeller
Cartidge passed the Germany
dynamic balance machine test

Imported bearing
Long life, low noise and stable

Reject welding,
complete head made
by CNC, more stable!

3 way spray,3 air
Head and cartridge with
anti-retraction structure

Professional LED light source,
superior brightness,Longevity
10000 hours.

Unique electric generator,
just a little air can generate
sufficient power;



This product has been applied layout-design
patent.counterfeit must be investigated.
NO.ZL201220717131.8 ZL201210563855.6

Technical:

High speed handpiece:	3 way spray 3 air
Chuck type:	push button
Bulb:	LED light
Noise	≤70dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥ 300,000rpm

Model: CX207-F

High-speed air turbine handpieces



Order code:**210-F+-SP**



Standard head

Order code: **H15-SP4**
H15-SP2



Order code:**210-F+-TP**



Torque head

Order code: **H15-TP4**
H15-TP2



Order code:**210-F-MP**



Miniature head Quick coupling

Handpiece order code: **H05-MPQ**
Complete order code: **H05-MPQ4**
H05-MPQ2



Order code:**210-F+-SP**



Standard head Quick coupling

Handpiece order code: **H15-SPQ**
Complete order code: **H15-SPQ4**
H15-SPQ2



Order code:**210-F+-TP**



Torque head Quick coupling

Handpiece order code: **H15-TPQ**
Complete order code: **H15-TPQ4**
H15-TPQ2

78
LED handpiece
with generator

LED HANDPIECE WITH GENERATOR

Model: CX207-F

Caries
detection
function

Self-luminous

Balanced impeller
Cartidge passed the Germany
dynamic balance machine test

Imported bearing
Long life, low noise and stable

Reject welding,
complete head made
by CNC, more stable!

3 way spray,3 air
Head and cartridge with
anti-retraction structure

LED for caries detection
Treat while checking
Faster and better

Unique electric generator,
just a little air can generate
sufficient power;



This product has been applied layout-design
patent,counterfeit must be investigated.
NO.ZL201220717131.8 ZL201210563855.6

Technical:

High speed handpiece:	3 way spray 3 air
Chuck type:	push button
Bulb:	LED light
Noise	≤70dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole)
	0.20-0.22MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥ 300,000rpm

Model: CX207-F

High-speed air turbine handpieces



Order code:210-F+-SP



Standard head

Order code: **H65-SP4**
H65-SP2



Order code:210-F+-TP



Torque head

Order code: **H65-TP4**
H65-TP2



Order code:210-F+-SP



Standard head Quick coupling

Handpiece order code: **H65-SPQ**
Complete order code: **H65-SPQ4**
H65-SPQ2



Order code:210-F+-TP



Torque head Quick coupling

Handpiece order code: **H65-TPQ**
Complete order code: **H65-TPQ4**
H65-TPQ2

79
LED handpiece
with generator

80
LED handpiece
with generator

OPTICAL HANDPIECE

Optical fiber illumination

Model: CX207-G

New upgrade, new integrated head



Fiber optic

Reject welding,
complete head made
by CNC, more stable!

Imported bearing
Long life, low noise and stable

Balanced impeller
Cartidge passed the Germany
dynamic balance machine test

3 way spray, 3 air
clean head system
Head and cartridge with
anti-dust structure

Optical fiber illumination

Graceful body, simple, portable
easy clean, easy control

Technical:

High speed handpiece:	3 way spray 3 air
Chuck type:	push button
Bulb:	LED light
Noise	≤70dB
Air pressure:	0.25-0.28 MPA(6hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥300,000rpm
Bulb power:	1w
Work voltage:	2.9v-3.2v
Work current:	160mA.

Model: CX207-G

High-speed air turbine handpieces



(Fit into **NSK**
Machlite/Phatelus*)

Standard head Handpiece Order code: **H16-NSPQ**
Complete order code: **H16-NSPQ6**

Torque head Handpiece Order code: **H16-NTPQ**
Complete order code: **H16-NTPQ6**



(Fit into **W&H**
Roto Quick*)

Standard head Handpiece Order code: **H16-WSPQ**
Complete order code: **H16-WSPQ6**

Torque head Handpiece Order code: **H16-WTPQ**
Complete order code: **H16-WTPQ6**



(Fit into **Sirona***R/F)

Standard head Handpiece Order code: **H16-SSPQ**
Complete order code: **H16-SSPQ6**

Torque head Handpiece Order code: **H16-STPQ**
Complete order code: **H16-STPQ6**



(Fit into **KaVo**
MULTIfelix LUX)

Standard head Handpiece Order code: **H16-KSPQ**
Complete order code: **H16-KSPQ6**

Torque head Handpiece Order code: **H16-KTPQ**
Complete order code: **H16-KTPQ6**



Order code: **210-G+-TP**
For torque head



Order code: **210-G+-SP**
For standard head



Order code: **229-GN**
(Fit into **NSK**
Machlite/Phatelus*)



Order code: **229-GW**
(Fit into **W&H**
Roto Quick*)



Order code: **229-GS**
(Fit into **Sirona***R/F)



Order code: **229-GK**
(Fit into **KaVo**
MULTIfelix LUX)

1 WAY SPRAY
HANDPIECE

Model: CX207-B

New upgrade, new integrated head

1 way spray, 4 air spray
Head and cartridge with
anti-dust structure

Reject welding,
complete head made by CNC,
more stable!

Good cooling system,
soft push button

Balanced impeller
Cartridge passed the Germany
dynamic balance machine test

Imported bearing
Long life, low noise and stable

Graceful body, simple,
easy clean, easy control

Wear-resisting, efficient operation





Technical:

High speed handpiece:	1 way spray, 4 air spray
Chuck type:	push button
Noise	≤68dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22 MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥ 300,000rpm

Model: CX207-B

High-speed air turbine handpieces



Order code: 210-B+-S



Order code: 210-B+-T



Order code: 210-B+-SP



Order code: 210-B+-TP



Standard head Order code: H13-S4
H13-S2



Torque head Order code: H13-T4
H13-T2



Standard head Order code: H13-SP4
H13-SP2



Torque head Order code: H13-TP4
H13-TP2

83
1 way spray
handpiece

84
1 way spray
handpiece

3 WAY SPRAY
HANDPIECE

Model: CX207-W

New upgrade, new integrated head

Cartridge passed the Germany dynamic balance machine test

3 way spray, 3 air spray
Head and cartridge with anti-dust structure

Good cooling system,
soft push button

Wear-resisting, efficient operation


Graceful body, simple,
easy clean, easy control

Technical:

High speed handpiece:	3 way spray, 3 air spray
Chuck type:	push button
Noise	≤68dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22 MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥300,000rpm

Model: CX207-W

High-speed air turbine handpieces



Order code: 210-W+-SP



Standard head Order code: H17-SP4
H17-SP2



Order code: 210-W+-TP



Torque head Order code: H17-TP4
H17-TP2



Order code: 210-W+-SP



Standard head Quick coupling Order code: H17-SPQ4
H17-SPQ2



Order code: 210-W+-TP



Torque head Quick coupling Order code: H17-TPQ4
H17-TPQ2

COLOR
HANDPIECE

Model: CX207-C

New upgrade, new integrated head

Ceramic bearing
Cartridge passed the Germany
dynamic balance machine test

Wear-resisting, efficient operation

Good cooling system,
soft push button



Aluminum alloys body, light,
colorful, easy clean, easy control

Technical:

High speed handpiece:	1 way spray, 4 air spray
Chuck type:	push button
Noise	≤68dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥300,000rpm
color:	red yellow purple blue black

Model: CX207-C

High-speed air turbine handpieces



Standard head
Red

Order code: H14-1SP4
H14-1SP2



Torque head
Red

Order code: H14-1TP4
H14-1TP2



Standard head
Yellow

Order code: H14-2SP4
H14-2SP2



Torque head
Yellow

Order code: H14-2TP4
H14-2TP2



Standard head
Purple

Order code: H14-3SP4
H14-3SP2



Torque head
Purple

Order code: H14-3TP4
H14-3TP2



Standard head
Blue

Order code: H14-4SP4
H14-4SP2



Torque head
Blue

Order code: H14-4TP4
H14-4TP2



Standard head
Black

Order code: H14-5SP4
H14-5SP2



Torque head
Black

Order code: H14-5TP4
H14-5TP2



Order code: 210-C+-SP
For standard head



Order code: 210-C+-TP
For torque head

Model: CX207-A

High-speed air turbine handpieces

New upgrade, new integrated head




Single way spray
more three air spray
Head and cartridge with
anti-dust structure

Cartridge passed the Germany
dynamic balance machine test

The special clean head designed
to avoid the entry of oral fluids and
other impurity into handpiece head.

Graceful body, simple,
easy clean, easy control



Order code: 210-A-SP



Standard head

Order code: H02-SP4
H02-SP2



Order code: 210-A-TP



Torque head

Order code: H02-TP4
H02-TP2

Model: CX207

High-speed air turbine handpieces

Power increased by more than 30%, cutting faster!

Bearing is imported
from Japan

Push button design
Cartridge is easier to replace

The neck is small and the treatment
area is maximized.






Technical:

High speed handpiece:	1 way spray
Chuck type:	push button
Noise	≤65dB
Air pressure at the back end of HP:	0.22-0.24MPa(4hole)
Effective range of overload pneumatic protection :	0.3MPa~0.5MPa
Rotation speed:	≥350,000rpm

Bearing is imported from Japan

Power increased by more than 30%,
cutting faster

Overload protection can prevent the working
pressure being too high and prolong the
service life of the handpiece



Order code: 210-1-SP



Standard head Order code: H 11- SP4

OPTICAL HANDPIECE

Optical fiber illumination

Model: CX207-G

High-speed air turbine handpieces



Order code: **210-GKD-SP**



Order code: **229-GK1**
fit into KaVo® MULTfelx® LUX

Standard head 45°

Handpiece Order code: **H06-KD1SPQ**
Complete order code: **H06-KD1SPQ6**



LED HANDPIECE WITH GENERATOR

Model: CX207-F

High-speed air turbine handpieces



Order code: **210-FD-SP**




Standard head 45°

Order code: **H05-DSP4**
H05-DSP2

Model: CX207

High-speed air turbine handpieces



- High speed handpiece 45 degree angle for impacted tooth.
- Water cooling , no air come out from the head , avoiding hematoma during treatment
- Small trauma, short time, accurate force
- High torque, stable speed
- Minimally invasive extraction of impacted teeth is faster and safer



Technical:

High speed handpiece:	1 way spray
Chuck type:	push button/latch
Noise	≤68dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥ 300,000rpm




Order code: **210-D-SP**



Standard head 45° Order code: **H01-D1SP4**
H01-D1SP2



Push button or wrench choose



Technical:

High speed handpiece:	1 way spray
Chuck type:	push button/latch
Noise	≤68dB
Air pressure at the back end of HP:	0.25-0.28 MPA(4hole) 0.20-0.22MPA(2hole)
Bur applicable:	Ø1.59-Ø1.6mmx21-23mm
Rotation speed:	≥ 300,000rpm



Order code:210-S



Order code:210-SP




Order code:210-TP



Standard head Order code: **H01-S4**
H01-S2



Standard head Order code: **H01-SP4**
H01-SP2



Torque head Order code: **H01-TP4**
H01-TP2

SPEED INCREASED
CONTRA ANGLE

Model: CX235C7

1:5 Speed increased Contra Angle

New upgrade, Better structure
Better performance

1:5 Speed increased
60° contra angle


1:5 speed increase

Technical:

Chuck type:	push button
Spray:	Inner channel
Bur applicable:	Ø1.6mm
Rotation speed:	about 200,000rpm

Model: CX235C7

Geared Angle Handpieces



Order code: C7-5
Optic fiber




External water
Order code: C7-5S
Optic fiber



Order code: C7-6



V shape water spray structure
More stable and better cooling




More concentrated light
Increase the illumination
by 100%



Add sanitary hole
Reduce the risk of dirt
entering



Better waterproof structure
To better prevent water from
entering the head



Improved shaft structure
Service life increased by
more than 90%

1:5 Speed increased
60° contra angle

SPEED INCREASED
CONTRA ANGLE

Model: CX235C7

1:5 Speed increased Contra Angle

New upgrade, smaller size
better handfeeling

Mini

97
Speed increased
contra angle




1.5 speed increase

Technical:

Chuck type:	push button
Spray:	Inner channel
Bur applicable:	Ø1.6mm
Rotation speed:	about 200,000rpm

Model: CX235C7

Geared Angle Handpieces



Reducing size and better angle

Improved handpiece and head
angle, making it easier to reach
the back teeth




98
Speed increased
contra angle



Smaller head with 92° angle

Larger elevation angle, making it easier to reach
the back teeth



Order code: C7-4

SPEED INCREASED
CONTRA ANGLE

Model: CX235C7

Contra angle handpiece

1:4.2 speed increase

Technical:

Chuck type:	push button
Spray:	Inner channel
Bur applicable:	Ø1.6mm
Rotation speed:	about 168,000rpm

4 hole water spraying
Spraying water to prevent subcutaneous emphysema during surgery
Effectively and quickly cool the needle.



96 Speed increased
contra angle

Perfect accessibility of 45 degree angle fiber-optic handpiece provides the best adaptability and maximum performance for teeth
Easy access to the molar area
Minimize interference with front and adjacent teeth
45 degree inclined head for supine accessibility and maximum visibility



Order code: C7-3
Optic fiber



External water
Order code: C7-3S
Optic fiber

SPEED INCREASED
CONTRA ANGLE

Model: CX235C7

Geared Angle Handpieces

1:5 speed increase

Technical:

Chuck type:	push button
Spray:	Inner channel
Bur applicable:	Ø1.6mm
Rotation speed:	about 200,000rpm



Stale, quiet and large torque cartridge

Skillful workmanship, reliable in quality

Specially for the European style, comfortable feeling

100 Speed increased
contra angle



optic fiber



Order code: C7-1
Optic fiber



Order code: C7-2

REDUCTION
CONTRA ANGLE

Model: CX235C6

Contra angle handpiece

20:1 Contra Angle for Implant

Easily disassembled, more reassuring treatment

101
Contra angle
handpiece



20:1 reduction

Technical:	
Chuck type:	push button
Spray:	Inner channel and external
Bur applicable:	Ø2.35 _{+0.016} mm(ISO1797-1)
Max speed :	2,000RPM
	For implant
Internal cooling system: Kirschner / Meyer	

Model: CX235C6

Geared Angle Handpieces



Easy to disassemble without tool
C6-22 can be disassembled with a simple twist for easy internal cleaning.however, when attached to a motor, unique locking mechanism makes sure the handpiece will not disasseb

102
Contra angle
handpiece



Optic fiber Order code: C6-22



Order code: C6-23

EXTERNAL
LOW SPEED HANDPIECE
SET

Geared Angle Handpieces



External water spray,
Short motor,
Reducing size and weight
More comfortable to use

External Order code: **235-M SET-4** ☺



Reducing size and weight
More comfortable to use

Technical:

Chuck type:	push button
Noise	≤70dB
Air pressure:	0.30MPA
Bur applicable:	Ø2.35 _{-0.012} mm(ISO1797-1)
Rotation speed:	about20,000rpm

EXTERNAL
LOW SPEED HANDPIECE
SET

Geared Angle Handpieces



External Order code: **235-F SET-4** ☺
235-F SET-2 ☺




External Order code: **235-A SET-4** ☺
235-A SET-2 ☺

INNER CHANNEL
LOW SPEED HANDPIECE
SET

Geared Angle Handpieces



Optic fiber inner channel Order code: 235-C SET-6



Optical fiber



6 hole



360 degrees of freedom rotary interfa

Technical:

Low speed handpiece: 1 way spray	
Chuck type:	push button
Noise	≤70dB
Air pressure:	0.30MPA
Bur applicable:	Ø2.35 _{-0.010} mm(ISO1797-1)
Rotation speed:	about20,000rpm

INNER CHANNEL
LOW SPEED HANDPIECE
SET

Geared Angle Handpieces



inner channel Order code: 235-B SET-4
235-B SET-2



Generator in contra angle



inner channel LED Order code: 235-E SET-4
235-E SET-2

Technical:

Low speed handpiece: 1 way spray	
Chuck type:	push button
Noise	≤70dB
Air pressure:	0.30MPA
Bur applicable:	Ø2.35 _{-0.010} mm(ISO1797-1)
Rotation speed:	about20,000rpm

REDUCTION CONTRA ANGLE

4:1 Reduction

Low speed handpiece:1:1 direct drive
Chuck type: push button or wrench
Spray: external
Bur applicable: Ø2.35_{±0.016}mm
Rotation speed: MAX10,000rpm



4:1 Order code: **CB-3**



Order code: **CH-1**



Order code: **CH-4**



Order code: **CH-11**



Up and down movement



Order code: **C3-1**




Order code: **C3-4**



Order code: **C3-11**

6:1 Reduction


Low speed handpiece: **6:1 Reduction**
Chuck type: push button
Max speed: 1200rpm
Max torque: 5.2N.cm



Order code: **C5-1M**

10:1 Reduction

Low speed handpiece:10:1 Reduction
Chuck type: push button
Spray: external
Bur applicable: engine use file/hand use file
Rotation speed: MAX4,000rpm
Rotation style: 90° reciprocating rotate head



Fit to hand use files



10:1 Order code: **CB-5**



Order code: **CH-12**



hand use files




Order code: **C5-12**

Geared Angle Handpieces

REDUCTION CONTRA ANGLE

64:1 Reduction

Low speed handpiece:64:1 Reduction
Chuck type: push button or wrench
Spray: external
Bur applicable: Ø2.35_{±0.016}mm
Rotation speed: MAX 625 rpm



64:1 Order code: **CB-8**



Order code: **CH-2**



Order code: **CH-4**



Order code: **C8-2**



Order code: **C8-4**

1:1 DIRECT DRIVE CONTRA ANGLE

1:1 Direct drive

Low speed handpiece:1:1 direct drive
Chuck type: push button or wrench
Spray: external
Bur applicable: Ø2.35_{±0.016}mm
Rotation speed: MAX40,000rpm



1:1 Order code: **CB-1**



Order code: **CH-2**



Order code: **CH-4**



Order code: **C1-2**



Order code: **C1-4**



1:1 Order code: **CB-7**



Order code: **CH-17**



Order code: **C-1F**

Geared Angle Handpieces

1:1 DIRECT DRIVE
CONTRA ANGLE

New upgrade, new integrated head

1:1 Direct drive

Low speed handpiece: 1:1 direct drive
Chuck type: push button
Spray: Inner channel
Bur applicable: $\varnothing 2.35_{-0.01}^{+0.02}$ mm
Rotation speed: MAX40,000rpm



Model: CX235-1C
Order code: C-1C
LED LIGHT



Model: CX235-1E
Order code: C-1E
GENERATOR



Model: CX235-1B
Order code: C-1B



Model: CX235C1
Order code: C1-1
Rotation speed: about 20,000rpm

AIR MOTOR

Spray: Inner channel or External
Hole type: 6 hole or 4 hole or 2 hole
Rotation speed: 20,000rpm



INNER CHANNEL

Model: CX235-3C
Order code: M-3C
Optic fiber



Model: CX235-3B
Order code: M-3B-4
M-3B-2



EXTERNAL

Model: CX235-3F
Order code: M-3F-4
M-3F-2



External

Model: CX235-3F
Order code: M-3M4

Geared Angle Handpieces

INNER CHANNEL/EXTERNAL
STRAIGHT HANDPIECE

1:1 Direct drive

Low speed handpiece: 1:1 direct drive
Spray: Inner channel or External
Bur applicable: For surgical $\varnothing 2.35_{-0.01}^{+0.02}$ mm
Rotation speed: MAX40,000rpm



Order code: S-2S
surgical operation 20°
straight head



Order code: S-2S2
surgical operation
straight head



INNER CHANNEL

Model: CX235-2C
Order code: S-2C
Optic fiber



Model: CX235-2B
Order code: S-2B



Model: CX235-2F
Order code: S-2F



EXTERNAL

Model: CX235-2A
Order code: S-2A






Model: CX235-2
Order code: S-2

ECONOMICAL
COMPLETE KIT

1:1 Direct drive


Economical complete kit



Order code: **H01-1 KIT-SP4** 
H01-1 KIT-SP2 



STUDENT HANDPIECE KIT

- 1PCS Anti-retraction high speed standard push handpiece
- 1PCS contra angle
- 1PCS straight head
- 1PCS low speed air motor
- 1PCS push button cartridge <Free>



111 Economical complete kit





Order code: **H01-2 KIT-S4** 
H01-2 KIT-S2 

HANDPIECE KIT


- 2PCS high speed standard handpiece
- 1PCS contra angle
- 1PCS straight head
- 1PCS low speed air motor
- 1PCS standard cartridge <Free>





Order code: **H01-2 KIT-SP4** 
H01-2 KIT-SP2 

HANDPIECE KIT

- 2PCS standard push handpiece
- 1PCS contra angle
- 1PCS straight head
- 1PCS low speed air motor
- 1PCS push button cartridge <Free>



Order code: **H01-1 KIT-S4** 
H01-1 KIT-S2 

HANDPIECE KIT

- 1PCS high speed standard handpiece
- 1PCS contra angle
- 1PCS straight head
- 1PCS low speed air motor
- 1PCS standard cartridge <Free>



112 Handpiece accessories

Handpiece Accessories
Manufacturing all kinds of spare parts for handpiece



OEM, ODM Service

Provide OEM,ODM Service

COXO has improved the design innovation system and a strong research and development capabilities, we can control the whole manufacture process according to the customers' requirements

HIGH FREQUENC AIR SCALER

MOQ:300PCS

Oem, odm service



Myonic bearing

LED Blue
3 way spray,3air
Head and cartridge with
anti-dust structure

Reject welding,
complete head made
by CNC, more stable!

Cartidge passed the Germany
dynamic balance machine test

Fit for **standard,Torque head**
and 45° degree handles,meet
with the dentist use perfectly

Unique quick coupling
with **the generator**

Separated generator with the handle,
free sterilization,
make the handpiece life longer



Order code:
210-F+-SP



Order code:
210-F+-TP



Order code:
210-FD-SP



Standard head
Order code:**H25-SP6**



Torque head
Order code:**H25-TP6**



Standard head 45°
Order code:**H25-DSP6**



Order code:**229-F4**
229-F2



Order code:**H25-3 KIT-4**
H25-3 KIT-2

LED handpiece with generator kit
• 1 PCS quick coupling for LED handpiece with generator(2/4 hole)
• 2 PCS standard push optical handpiece with generator
• 1 PCS torque push optical handpiece with generator



Order code:**H25-5 KIT-4**
H25-5 KIT-2

LED handpiece with generator kit
1 PCS quick coupling for LED handpiece with generator(2/4 hole)
2 PCS standard push optical handpiece with generator
1 PCS torque push optical handpiece with generator
1 PCS contra-angle block teetn standard push optical handpiece

CUSTOM PRODUCT
HANDPIECE

MOQ:500PCS



LED HANDPIECE WITH
GENERATOR

3-way spray,3-way air



(Fit into KaVo)



Standard head



Complete order code: **H35-SPQ**

Handpiece order code:**H35-SPQ4**




Torque head



Complete order code: **H35-TPQ**

Handpiece order code:**H35-TPQ4**

MOQ:300PCS



Order code:**K07**

Torque

(For KaVo®630PB)



Torque head

Order code: **H01-K1TP4**

H01-K1TP2



Order code:**W08**

Torque

(For W&H®TA-97)



Standard head

Order code: **H01-W1SP4**

H01-W1SP2



Order code:**W05**

Torque

(For W&H®Synea TA-98)



Torque head

Order code: **H01-W1TP4**

H01-W1TP2

MOQ:500PCS



(For W&H C-SWSABT/Trend)



LED HANDPIECE WITH
GENERATOR



Torque head

Order code: **H05-W2TP4**

H05-W2TP2



Torque head

Order code: **H01-W2TP4**

H01-W2TP2

CLEAN HEAD SYSTEM

MOQ:300PCS



LED HANDPIECE WITH
GENERATOR



Standard head

Order code: **H15-8SP4**

H15-8SP2



Torque head

Order code: **H15-8TP4**

H15-8TP2

LED HANDPIECE WITH
GENERATOR

MOQ:300PCS



LED HANDPIECE

Torque head Order code: **H26-GL**

CUSTOM PRODUCT
HANDPIECE

MOQ:500PCS

3 way spray 3 air spray handpiece

- A1-SP
Stan push head
- A1-TP
Max push head
- A1-45SP
Stan push head 45°

Clean head system handpiece

- A2-S
Stan head
- A2-T
Max head
- A2-SP
Stan push head
- A2-TP
Max push head

Three way spray handpiece


- A3-SP
Stan push head
- A3-TP
Max push head

LED handpiece with generator

- A4-SP
Stan push head
- A4-TP
Max push head
- A4-45SP
Stan push head 45°

Fiber optical handpiece

- A5-SP
Stan push head
- A5-TP
Max push head



B1
body and pattern



B2
body and pattern



B3
body and pattern



B4
body and pattern



B5
body and pattern



C1
Coupling




C2
Quick coupling(for NSK)



C3
Quick coupling(for Sinol)




C4
Quick coupling of fiber optical
(for NSK)




C5
Quick coupling of fiber optical
(for W&H)



C6
Quick coupling of fiber optical
(for Sirona)



C7
Quick coupling of fiber optical
(for KaVo)



C8
Quick coupling
(for KaVo)



2 hole




4 hole



6 hole

COXO could Provide OEM,ODM service,and COXO have some handpiece style and pattern provide customers choice, Choose as shown in figure:

customized order



Dental Operating Lights
Professional dental products manufacturer



Head Office: Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE

Phone: +971 4 569 3474 Mobile: +971 50 912 4577 Email: admin@skyline-med.com

