

Clinical Cabinets and Accessories



"BEYOND YOUR
EXPECTATIONS..."

Anemon Series





PRD. CODE AN00

CLINIC CABINET WITH FOOT CONTROL



Dimensions

(W x L x H) 500 x 2000 x 860. (± 50 mm)

Top: Corian

Sink: Corian

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Glass

Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Adjustable chrome legs

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE AN04

CLINIC CABINET



Dimensions

(W x L x H) 500 x 2200 x 860. (± 50 mm)

Top: Glass

Sink: Glass

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Corian

Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Adjustable chrome legs

15 mm Tempered Glass top

Soft-closing sliding system (Blum)

Iron phosphat coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE AN05

WALL MOUNT CLINIC CABINET



Dimensions

(W x L x H) 500 x 2100 x 860. (± 50 mm)

Top: Corian

Sink: Corian

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Glass

Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Adjustable chrome legs

15 mm Tempered Glass top

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE AN06

WALL MOUNT CLINIC CABINET



Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Adjustable chrome legs

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust

Dimensions

(W x L x H) 500 x 2000 x 860. (± 50 mm)

Top: Corian

Sink: Corian

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Glass



WWW.SKYLINE-MED.COM



PRD. CODE AN07

WALL MOUNT CLINIC CABINET



Specifications

- 0.90-1 mm metal sheet
- Electrostatic powder coating
- Various ral color choice
- Soft-closing sliding system (Blum)
- Iron phosphate coating against the rust

Dimensions

(W x L x H) 500 x 1500 x 860. (±50 mm)

(W x L x H) 500 x 500 x 860. (±50 mm)

Top: Corian

Sink: Corian

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Glass



"BEYOND YOUR
EXPECTATIONS..."

Azurit Series





PRD. CODE AZ01

CLINIC CABINET WITH MOBILE



Specifications

- 0.90-1 mm metal sheet
- Electrostatic powder coating
- Various ral color choice
- Adjustable chrome legs
- 12-24 mm Corian top
- Soft-closing sliding system (Blum)
- Iron phosphate coating against the rust

Dimensions

(W x L x H) 500 x 2100 x 860. (± 50 mm)

Top: Corian

Sink: Corian

Options

- Antibacterial plastic drawer organizers (ABS)
- Top Choices: Laminate / Glass
- Foot control for mixer (Optional)



WWW.SKYLINE-MED.COM



PRD. CODE AZ02

CLINIC CABINET



Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Adjustable chrome legs

12-24 mm Corian top

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust

Dimensions

(W x L x H) 500 x 2000 x 860. (± 50 mm)

Top: Corian

Sink: Corian


Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Laminate / Glass

Foot control for mixer (Optional)





PRD. CODE AZ05

CLINIC CABINET WITH PICTURE



Specifications

- 0.90-1 mm metal sheet
- Electrostatic powder coating
- Various ral color choice
- Various picture / logo painting
- Adjustable chrome legs
- 15 mm glass top
- Soft-closing sliding system (Blum)
- Iron phosphat coating against the rust

Dimensions


(W x L x H) 500 x 2000 x 860. (±50 mm)

Top: Glass

Sink: Glass

Options

- Antibacterial plastic drawer organizers (ABS)
- Top Choices: Laminate / Corian



WWW.SKYLINE-MED.COM

"BEYOND YOUR
EXPECTATIONS..."

Kuvars Series





PRD. CODE KV01

CLINIC CABINET



Dimensions

(W x L x H) 500 x 2100 x 730. (± 50 mm)

Top: Glass

Sink: Glass

Options

Antibacterial plastic drawer organizers (ABS)

Top Choices: Corian

Specifications

0.90-1 mm metal sheet

Electrostatic powder coating

Various ral color choice

15 mm tempered glass top

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE KVM

CLINIC CABINET WITH MOBILE



Specifications

1-4 mm metal sheet

Electrostatic powder coating

Various ral color choice

Silent castors

Corian organizer

15 mm glass top

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust

Dimensions

(W x L x H) 500 x 500 x 830. (± 50 mm)

Top: Corian

Options

Antibacterial plastic drawer organizers (ABS)



WWW.SKYLINE-MED.COM



"BEYOND YOUR
EXPECTATIONS..."

Mobile Cabinets





PRD. CODE MB02

MOBILE CABINET



Dimensions

(W x L x H) 500 x 500 x 830. (± 50 mm)

Top: Corian

Options

Antibacterial plastic drawer organizers (ABS)

Specifications

1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Silent castors

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE MB03

MOBILE CABINET



Dimensions

(W x L x H) 500 x 500 x 830 (± 50 mm)

Top: Glass

Options

Antibacterial plastic drawer organizers (ABS)

Specifications

1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Silent castors

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust



WWW.SKYLINE-MED.COM



PRD. CODE MB04

MOBILE CABINET



Dimensions

(W x L x H) 500 x 500 x 830 (± 50 mm)

Top: Laminate

Options

Antibacterial plastic drawer organizers (ABS)

Specifications

1 mm metal sheet

Electrostatic powder coating

Various ral color choice

Side instrument shelf

Silent castors

Soft-closing sliding system (Blum)

Iron phosphate coating against the rust




WWW.SKYLINE-MED.COM

"BEYOND YOUR
EXPECTATIONS..."

Lab.
Products





PRD. CODE PH01

DENTAL SIMULATION UNIT WITH PHANTOM

Specifications

High and low speed handpiece, suction and three-way syringe

It is suitable central air and water system

Unit type lighting is used as the light source

The body is electrostatically painted sheet metal

Laminate top

Dimensions

(W x L x H) 900 x 600 x 900 (h)





PRD. CODE GPW01

DENTAL WORKBENCH SINGLE

Specifications

80 kg. weight.

220 V / AC 50 Hz. / Max. 1250 Watt.

Produced of CNC manufactured, 1 mm D.K.P. steel Sheet.

Electrostatic powder coating.

Various RAL color range.

Gas tap

30mm laminate post-formed work with 3mm PVC side cover.

Air gun.

Holes for micromotor handpiece on worktop.

Dimensions

(W x L x H) 1050 x 650x 850mm (h)





PRD. CODE GPW02

DENTAL WORKBENCH DOUBLE

Specifications

160 kg. weight

220 V / AC 50 Hz. / Max. 1250 Watt

Produced of CNC manufactured, 1 mm D.K.P. steel Sheet

Electrostatic powder coating

Various RAL color range

Gas tap

30mm laminate post-formed work with 3mm PVC side cover

Air gun

Holes for micromotor handpiece on worktop

Dimensions

(W x L x H) 2050 x 650x 850mm (h)



Clinical Cabinets and Accessories



2020 COLLECTION

Clinical Cabinets and Accessories



Head Office: Office 509, Churchill Executive Tower,
Al Amal St., Business Bay, Dubai,UAE

Phone: +971 4 569 3474

Mobile: +971 50 912 4577

Email: admin@skyline-med.com

