

GENI

Dental Anesthesia Booster


Beautifully design
Comfortable to use



3 Adjustable Injection Speeds

- High-precision stable injection
- Doctors can choose different speed modes according to habits

Mode	Time to inject 1 cartridge of anesthetic
Low	About 285s
Mid	About 130s
High	About 89s



Soft touch control

- Induction injection, More convenient for doctors to use
- One-button switching mode, adjustable speed during the injection process

Music system

- Set up speed mode voice prompt, injection status voice prompt
- Comes with soothing music, which can relieve the tension of patients

