



Powder Jet Handpiece

Simple and Efficient Plaque Control Solution

ProJet

Simple and Efficient Plaque Control Solution

The new VOLDENT ProJet 2 meets all the needs of daily prophylaxis treatment. It's suitable for working in subgingival or supragingival regions, gaps, bridges, implants and crowns. It's feather-light, perfectly balanced and ergonomically designed for comfortable working.



ProJet 2

- 3 sandblasting modes adjustable
- Equipped with 2 sandblasting handpiece PJ-02 (supragingival) & PJ-01 (periodontal)
- KAVO MULTiflex / MILTiflex LED quick connector



ProJet 1P

- Equipped with PJ-01 supragingival sandblasting handpiece
- KAVO MULTiflex / MILTiflex LED quick connector



ProJet 1

- Equipped with PJ-01 supragingival sandblasting handpiece
- 4-holes connector



ProJet

Simple and Efficient Plaque Control Solution



Unique Periodontal nozzle
The precise little cannula for periodontal pockets

The narrow plastic periodontal sandblasting nozzle for focused gentle periodontal treatment as well as cleaning gaps, bridges, implants and crowns.

0.7mm powder outlet
Exceptional cleaning and concentrated jet

Less powder mist thanks to the highly precise power jet. The narrow shape of the cannula provides concentrated sandblasting during supragingival cleaning.

Enhanced visibility
Thin and long nozzle, pre-mixed powder and water

Work comfortable even in difficult-to-reach places, thanks to the thin and long nozzle which is 360° rotatable. The powder and water are pre-mixed before spraying, minimizes blowing powder.

Stay relaxed, handle with easy

A super light, perfectly balanced handpiece with individually adjustable lengths for comfortable handling and stress-free treatment.

3 power levels of adjustment

Three power levels adjustment, with the first level for powerless cleaning of the teeth to remove powder residue. Level 2 and 3 are for individually adjustable cleaning performances depending on the patient and indication.

3-Section fuselage design

Segmented handpiece minimizes the risk of clogging. This is the easiest sandblasting handpiece to unclog and clean.

Simple filling and safety storage

A second powder container for interruption-free treatment, the powder container can be filled in advance and sealed with a rubber seal which shows the indication.

Easy quick coupling

Can fit with both Kavo Multiflex LED and no-LED quick coupling.





Head Office: Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE
Phone: +971 4 569 3474 **Mobile:** +971 50 912 4577 **Email:** admin@skyline-med.com

