


# S-AOPT

*Brighter than a shooting star, Shine no matter where you are.*

**NINE** Wavelength  
Filters Optional

*Utilize OPT And MSP Technology*





# *With the S-APOT, The Skin Rejuvenated is Easy.*



Acne



Epidermis pigmentation



Epidermis pigmentation,  
Rosacea and flushing



Pigmentation,  
Vascularisation and skin care.



Dermis larger vessels,  
Skin care,  
Pigmentation and vascularization



Larger vessels in the face,  
General ipl skin care,  
Pigmentation and vascularization



Oil control,  
Fine lines,  
Skin texture improvement



Hair removal




Deep anti-aging




# Principle

*The longer the wavelength, the deeper the penetration*



# OPT and MSP Technology


## OPT - Optimal Pulse Technology



### The OPT advantages:

- High peak power, shorter pulses - optimal for IPL skin treatments using photorejuvenation and the treatment of benign pigmented lesions.
- Homogeneous power over long pulse durations.

## MSP - Multiple-Sequential Pulsing



### MSP Technology:

- Multiple-Sequential Pulsing allows for cooling the skin in between pulses, reducing the chance of skin damage during treatments.



# Cooling Handle

-4-20°C Cooling ❄️



# **Spot Size Adaptors**

## ***Quick-Change Spot Size Adaptors***

Quickly and easily change spot sizes with magnetic gold pieces that instantly adjust for comfortable treatment of both small and highly contoured areas.





# **Advantages**




***One versatile applicator to treat over 20 skin conditions and hair removal.***



**Treats Acne, Rosacea,  
Pigmented and Vascular  
Lesions.  
Also reduce age spots  
through photorejuvenation.**



**Optimal Pulse Technology  
(OPT) controls pulse shape  
and length.**



**Multiple Sequential Pulsing  
(MSP) enables a safer  
procedure while cooling  
between pulses.**



**9 Filters, changeable in  
seconds.**



**-4-20°C for patient  
comfort through  
continuous contact cooling.**

# UI Intertace



Pulse Duation 5-100ms

Interval Time Between  
Two Shots 1-5s

Sub Pulse

Auto Identified Filters

Energy 5-25j/cm<sup>2</sup>

Temperature -4°C-20°C  
Adjustable



# Parameter

## Specification of S-AOPT

### WAVE LENGTH

420-1200nm

### FILTERS

420nm;530nm; 515 nm, 560 nm; 590 nm; 615 nm;  
640 nm; 695 nm; 755 nm;

### FLUENCE

5-25J/cm<sup>2</sup>

### INPUT POWER

3000W

### POWER SOURCE

Pure Sapphire

### LAMP

UK Imported lamp

### SPOT SIZE

15\*50mm

### CONTACT COOLING TEMPERATURE

-4°C-20°C adjustable

### PACKAGE SIZE

65\*55\*121cm

### COOLING SYSTEM

Built in water cooling +semi conductor cooling +air cooling

### WEIGHT

87kg

### LCD SCREEN

10.4 inch true color touch screen control



# Shrink pores



Before

After Immediate effect

•We're confident in our results, so you can be too



# Rosacea



**Before**

**After 2 treatments**

•We're confident in our results, so you can be too

# Freckle



Before



After 2 weeks

•We're confident in our results, so you can be too



# **Vascular**



**Before**




**After 4 weeks**

•We're confident in our results, so you can be too

# Acne inflammation



Before




After 2 treatments

•We're confident in our results, so you can be too




# lass lip



Before

After Immediate effect

•We're confident in our results, so you can be too



# SKYLINE

**MEDICAL REQUISITES TRADING LLC**

**Head Office:** Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE  
**Phone:** +971 4 569 3474 **Mobile:** +971 50 912 4577 **Email:** admin@skyline-med.com

