

ES-20

Electrosurgery Unit


**Safe/Professional/Effective
Convenient/Reliable**

Fast and precise Incision and coagulation
of soft tissue

With 7 pieces high performance and
hyperfine electrodes

Excellent therapeutic control

Innovative conductive bracelet (Neutral Electrode)



Specifications:

Dimension of Main Unit: 23 * 26 * 11 cm

Weight of Main Unit: 1.5kg

Supply Voltage: 115VAC / 230VAC

Supply Frequency: 50Hz / 60Hz

Power Consumption: 170VA to 230VA

Max Power Output: 30W