

INDEEPLUS 448khz PRO

* CAPACITIVE AND RESISTIVE DIATHERMY with Vacuum *




THE REVOLUTION Technology At 448 kHz in



• Anti-wrinkle

• Anti-aging

• Rejuvenation



Fat reduction

Body contouring

Smoothness

Tightening

WHAT IS INDEEPLUS 448KHZ?

Indeeplus 448KHZ is a unique and patented radiofrequency, so called Proionic®radiofrequency, which offers pleasant and effective treatments for facial and body enhancements.

The latest innovation in aesthetic appliances INDEEPLUS 448KHZ is an incredibly versatile new solution that allows you to have a wide variety of treatments in a single device, at the 448kHz frequency.

Is indicated for **skin tightening, cellulite reduction, body contouring and pre/post-surgery applications.**

It is the most advanced solution on the market, because it puts the latest cutting-edge medical-aesthetic technology at your fingertips. A true alternative that allows you to go one step further, making your investment profitable in less time.





Increasing core body temperature up to 40-50°C



Expansion of artery and capillary



Increasing blood flow rate



Increase in Collagen

Stimulation of Immune system

Increase in Metabolism

Removal of toxins



Face lifting



Healthy




Obesity care



Reduce pain

RF diathermy



CAP (Capacitive Electronic Transfer)
3 to 5 cm depth deep



Resistive Diathermy treatment with the Vacuum system





The capacitive electrode focuses the therapy into muscle layers.

Aiming at the superficial tissues of 3-5cm of skin, the effect is firming and anti-aging, it helps to restore the elasticity of blood vessels, rebuild the microcirculation system, provide more nutrients for fibroblasts, strengthen collagen synthesis, improve cells the flow of reactive oxygen, thereby restoring cell viability and rebuild the reticular layer.



Resistive Diathermy treatment the Vacuum system

Vacuum pressure utilizes numerical control rhythm mode combined with specially designed negative air suction tip, according to the personalized skin condition of the human body, respectively, the skin's epidermis, dermis, and fat layers are kneaded and massaged at different depths to stimulate the body's cells effectively. The flow of body fluids increases cell movement, activates cells, promotes blood circulation and lymphatic circulation of tiny blood vessels, accelerates metabolism, and improves the internal environment of the skin



TREATMENT EFFECT


FACIAL TREATMENTS

Wrinkles · elasticity · eyelids · neck · cleavage



BODY TREATMENTS

Chest · sagging · fat · cellulite · stretch marks · adiposity





ADVANTAGE

1.The RF output power is more sufficient, the maximum power is 300W, and the energy rises faster.

2.Add the Vacuum system to your treatment Add to your Resistive Diathermy treatment the Vacuum system, which will help in the treatment by mobilizing fat and increasing drainage.

3.Dynamically adjust the depth of action of energy and automatically adjust the energy flow to ensure that the tissue can receive the most effective treatment, the treatment is more comprehensive, and the effect is better.

4.Perceive the type of treatment head to adjust the output waveform and energy to make the output more stable and controllable.









5.Real-time perception of body circuit impedance, thereby adjusting the output, making the experience more comfortable.

SPECIFICATIONS

Name	Indeepus RF 448khz
Frequency	448KHZ
Treatment handles	CAP y RES with vaccm
Output power	300W
Monopolar Capacitive Electrode:	φ40, 50, 60, 70mm
Monopolar resistive electrode	φ30, 50, 70, 80mm
Voltage	220V 50HZ o 110V 60HZ
Machine dimensions	320*370*820 mm
Net weight	15.5kg
Package	560*500*980(mm)



HANDLE



INTERFACE



FACE TREATMENT

For facial treatment, it helps regenerate the skin structure by increasing the synthesis of collagen and elastic fibers, so that unstructured tissues are hydrated and oxygenated, thereby improving its quality, thickening the skin and giving it greater resistance.




Reduce fine wrinkles
Improve eye bags and dark circles
Facial rejuvenation
Reduce double chin

One course contains 6 treatment times
Each treatment time around 40 minutes
Interval time 1-2 weeks



BODY TREATMENT

For body treatment, it can rebalance the cell exchange of ions and increase metabolism, thereby oxygenating tissues and helping cells restore their natural functions. As a result, the skin tissue regains the lost elasticity, the fat gradually disappears, the correct microcirculation is restored, and overall, the skin is rejuvenated.



Fat reduction

- **Body reshaping**
- **Cellulite treatment and partial skin reconstruction.**
- **Increases the level of collagen**
- **Tighten loose skin**
- **Muscle relaxation, pain relief.**



TREATMENT COURSE



How many sessions of INDEEPLUS treatment do we need?

Depending on the indications and areas of treatment, a range of 5-6 or 8-10 sessions may be necessary. The number of sessions will be adapted according to the medical prescription and individual patient need.

How frequently treatments can be undertaken?

Usually, there are 1 to 2 treatment sessions a week but there is no limitation. As the radiofrequency 448khz is completely harmless for the human body and its biological tissues. The frequency of the sessions will be adapted according to patient need.

A woman with long dark hair is lying on her back, looking upwards. A treatment device is applied to her abdomen, emitting a bright light. The background is a soft, out-of-focus grey.

WHAT ARE THE POSSIBLE SIDE EFFECTS PRODUCED BY INDEEPLUS 448KHZ?

There are no side effects as the Proionic® radiofrequency is innocuous for the human body and its biological tissues. However on occasion some very sensitive skin may experience some redness and tingling.

THE TREATMENT IS PAINFUL?


No, the treatment gradually and homogeneously increases temperature at selected areas in the body while respecting the cell physiology. It therefore produces a pleasant and non-invasive treatment without adverse effects.

WHAT ARE THE CONTRAINDICATIONS?

- Pacemaker or other electronic implants
- Pregnancy
- Non-intact skin (open wounds or burns)
- Thrombophlebitis

Health Beyond Health

INDEEPLUS





Head Office: Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE
Phone: +971 4 569 3474 **Mobile:** +971 50 912 4577 **Email:** admin@skyline-med.com

