




SKYLINE
MEDICAL REQUISITES TRADING LLC



INDEEPRO



Anti-wrinkle



Anti-aging



Rejuvenation

FACE

BODY



Fat reduction



Body contouring



Smoothness



Tightening



Increasing core body temperature up to 40-50°C



Expansion of artery and capillary



Increasing blood flow rate



Face lifting

Healthy

Obesity care

Reduce pain

Increase in
Collagen

Stimulation of
Immune system

Increase in
Metabolism

Removal of
toxins


RF DIATHERMY



CET



BIPOLAR



RET



CAP (Capacitive Electronic Transfer)
3 to 5 cm depth deep

Resistive Diathermy treatment with
the Vacuum system




RET



RET


RET working mode suitable for thicker, fatty and fibrotic tissues, good for big area treatment like body





CET

CET working mode suitable for superficial tissues, good for small area treatment like face, skin tightening and anti-aging effects.



HSP +



HSP Positive

Fixed area treatment without sliding treatment

HSP -





HSP Negative

Fixed area treatment without sliding treatment

ADVANTAGE

- 1.The RF output power is more sufficient, the maximum power is 300W, and the energy rises faster.
- 2.Add the Vacuum system to your treatment Add to your Resistive Diathermy treatment the Vacuum system, which will help in the treatment by mobilizing fat and increasing drainage.
- 3.Dynamically adjust the depth of action of energy and automatically adjust the energy flow to ensure that the tissue can receive the most effective treatment, the treatment is more comprehensive, and the effect is better.
- 4.Perceive the type of treatment head to adjust the output waveform and energy to make the output more stable and controllable.
- 5.Real-time perception of body circuit impedance, thereby adjusting the output, making the experience more comfortable.

UI INTERFACE



SPECIFICATIONS

NAME	INDEEPRO
Treatment handles	RET、 CET、 HSP
RET	R70、 ϕ 70mm R50、 ϕ 50mm R40、 ϕ 40mm R30、 ϕ 30mm
CET	C70、 ϕ 70mm C50、 ϕ 50mm C40、 ϕ 40mm C30、 ϕ 30mm
Frequency	448KHZ
HSP	143*86mm
Output power	500W
Voltage	220V 50HZ or 110V 60HZ
Machine dimensions	318*294*498mm
Net weight	14.6KG
Package	480*385*695mm

BEFORE & AFTER



Before

After



Before

After



Before

After



Before

After



Before

After



Before

After