P14D

Technical specifications

Intended use	Endodontics						
Indications for use	Retro preparation of cavities and isthmuses						
Description	Diamond-coated universal retro tip for preparing canals of anterior teeth						
Tip type	Smooth						
Allergenic factor	Tip made of su 12b Surgical stainle • nickel • chromiu						
Commercial reference	P14D : F00106						
Medical class	Class of medical device: IIa according to 93/42/EEC directive						
Adjustments	Newtron range	P5 Range	Suprasson P5 Booster		Prophy Max - P-Max - Implant Center 2 - Piezotome 2		
	7	7	5	3		•••	

Service life

As it is not possible to establish a maximum number of uses (that may be determined by many parameters such as duration of use, hardness of tissue, the force applied, wear), we recommend that medical devices being routinely used are replaced at least once a year. Replace the medical device if the tip does not vibrate at the expected frequency, the treatment is not progressing as normal and is taking longer or is at a standstill.

Associated documentation


Document title	References
Consulting electronic user instructions	J00007
Cleaning, disinfection and sterilisation instructions for tips	J02001
Apical Surgery clinical tips	J08061
Ultrasonic generator power settings table	J58000
General instructions relating to the complete range of standard ultrasonic dental tips	J02101




P14D

Sterilisation

Sterilise prior to each use



The tips must be cleaned, disinfected and sterilised prior to each use. Prior to each use, please refer to the cleaning, disinfection and sterilisation instructions for tips J02001.






The practitioner has to continually monitor the clinical site to immediately identify any heating risk.

Tip | P14D | J00401 | (98) | V4 | 05/2021 | NB17EN010D

SATELEC|A Company of ACTEON Group 17 av. Gustave Eiffel| ZI du Phare 33700 MERIGNAC|FRANCE Tel. +33 (0) 556 34 06 07|Fax. +33 (0) 556 34 92 92 satelec@acteongroup.com|www.acteongroup.com

Installation



Manufacturer identification



SATELEC

A Company of ACTEON Group 17, avenue Gustave Eiffel ZI du Phare 33700 MERIGNAC France Tel. +33 (0) 556.34.06.07 Fax. +33 (0) 556.34.92.92 E.mail : satelec@acteongroup.com www.acteongroup.com



Date of first CE marking

1998



