



Surf the canal with confidence

Wave • one[®]
GOLD



One complete solution

With WaveOne® Gold, DENTSPLY Tulsa Dental Specialties promotes a comprehensive treatment solution and optimizes predictable, successful outcomes



Keeping the WaveOne® simplicity with a golden host of benefits


wave • one®
GOLD

1

Reinforces patient safety

- Primary WaveOne® Gold file is 50% more resistant to cyclic fatigue than WaveOne Primary file
- Reduced screwing effect compared to standard rotary systems

Cyclic fatigue resistance




2

Covers a wider range of canal morphologies

- Enhanced file flexibility thanks to the Gold technology
- Extended size range (Small, Primary, Medium, Large)

Flexibility




3

Shortens the shaping time

- Takes the cutting efficiency to a higher level
- A single file per treatment also translates into time savings for canal shaping and irrigation

Shaping working time



Engineered with the latest reciprocating technology

DENTSPLY Tulsa Dental Specialties has teamed up with the best endodontic professionals to develop WaveOne® Gold

The design team speaks




Dr.C. Ruddle, Dr.J. Webber, Dr.W. Pertot, Dr.S. Kuttler

Shaping canals has never been easier.


Since its launch in 2011, DENTSPLY Tulsa Dental Specialties' WaveOne® has set a new standard in root canal shaping. All around the world general dentists and endodontists have voiced their appreciation of the simplicity WaveOne brings to their efforts to provide quality endodontic treatment to their patients.

We have been working on a new and exciting WaveOne project to bring you WaveOne Gold.

WaveOne Gold brings to the endodontic procedure a metallurgically advanced single-file technique for shaping canals. Further, we optimized the tip diameters, tapers, and the cross-section to produce a file that really improves safety, efficiency and flexibility when preparing canals. WaveOne Gold takes simplicity and treatment success to another level.




Courtesy of Dr.J. Webber



Watch the video on [YouTube](#) [f](#)

Simplicity that enables treatment focus

The majority of root canals can be shaped with a single WaveOne® Gold Primary file



Ready-to-use + Single use



Most dental practitioners have the twofold aim of reducing overall treatment costs while improving safety and efficiency. With this in mind, WaveOne Gold files are supplied in sterilized blister packs. Dentists will benefit from:

- Simplicity, thanks to the ready-to-use sterile packaging
- Increased safety margin, by always using files with maximum fatigue resistance potential
- Reduced risk of cross-contamination, which can occur when instruments are reused
- Maximized cutting efficiency

Product references

wave•one[®]
GOLD

Reciprocating Files



small reciprocating file (3 per pack)

WGSMA21 (20/.07) – 21 mm

WGSMA25 (20/.07) – 25 mm

WGSMA31 (20/.07) – 31 mm



primary reciprocating file (3 per pack)

WGPRIME21 (25/.07) – 21 mm

WGPRIME25 (25/.07) – 25 mm

WGPRIME31 (25/.07) – 31 mm



medium reciprocating file (3 per pack)

WGMEDIUM21 (35/.06) – 21 mm

WGMEDIUM25 (35/.06) – 25 mm

WGMEDIUM31 (35/.06) – 31 mm



large reciprocating file (3 per pack)

WGLARGE21 (45/.05) – 21 mm

WGLARGE25 (45/.05) – 25 mm

WGLARGE31 (45/.05) – 31 mm

Starter Kits



ProMark[®] & e3[®] Starter Kits:

7 Packs WaveOne[®] Gold Reciprocating Files
– any size

1 Motor (ProMark[®] or e3[®])

3 Packs ProGlider[®] Rotary Glide Path Files
(3 files per pack, any length)

8 Packs Pre-sterilized Lexicon[®] K-Files
(6 files per pack, any size)

1 Endo Access Kit (Autoclavable Holder with
8 Access Instruments)

1 Box ProLube[®] Root Canal Conditioner
(35 ampoules per box)

3 Boxes ProRinse[®] Endo Irrigation Needles
(5 needles per box)

1 Bottle QMix[®] 2in1 Irrigating Solution
(60 mL)

Obturator




GuttaCore[®] (5 per pack) (25 per pack)

yellow - GCP02026 GCP02023

red - GCP02526 GCP02523

blue - GCP03026 GCP03023

black - GCP04026 GCP04023




Obturator (6 per pack) (30 per pack)

small - WGOB6S WGOB30S

primary - WGOB6P WGOB30P

medium - WGOB6M WGOB30M

large - WGOB6L WGOB30L



Size Verifiers




Size verifiers (6 per pack)

small - WGSVS

primary - WGSVP

medium - WGSVM

large - WGSVL



Gutta-percha points



Gutta-percha points (60 per box)

small - WGGPCS

primary - WGGPCP

medium - WGGPCM

large - WGGPCL



Paper points



Sterile paper points (180 per box)

small - WGAPSS

primary - WGAPSP

medium - WGAPSM

large - WGAPSL



Equipment



ProMark[®] Motor - PMKEM1

e3[®] - E3KIT

GuttaCore[®] Oven - HTR-GC

WaveOne[®] Oven - HTR-WO



1-800-662-1202

For the latest information consult
www.TulsaDentalSpecialties.com

DENTSPLY
TULSA DENTAL
SPECIALTIES