

**17 IN 1**

# *Synthetic Instrument*



Lead-in



Hydration



Lifting



Anti-aging



Rejuvenation



Therapy



Improved



Recovery





Anti-ageing  
Firming



Deep Cleansing  
Nourish Skin



Therapy  
Recovery

Hydration  
Water Retention



Non-invasive  
Transdermal Absorpt



Brightens  
Rejuvenates




## Hydra Dermabrasion Machine

Deeply cleanses the skin system, cleansing blackheads, acne, and concealing insects; while replenishing moisture, nutrients; refines pores and increases skin clarity; improves skin texture; led therapy.

# Say Goodbye to Dull or Yellow Skin

Our hydration equipment reaches right into the dermis, replenishing and locking in water to restore the skin's normal water cycle.

## 1-Corneal Layer Dehydrated




Small molecules of skin internal water shortage, the skin's ability to protect the skin will be worse, the corresponding skin will have some sensitive, rough and other problems.

## 2-Dehydration of the Granular Layer



Functions such as protection against ultraviolet rays will be weakened or lost, then the skin's ability to resist melanin will also disappear, so that the skin will be prone to blemishes or make acne marks in the skin deepen in colour.

## 3-Basal Layer Dehydration



Cells in the basal layer are not able to grow normally because of the lack of water, so the metabolism of the skin will be slow, resulting in the deposition of pigments in the skin, which cannot be metabolised, resulting in yellowish skin or even blemishes and other skin problems.

## 4-Dermal Dehydration



The entire water cycle will have problems, which will lead to wrinkles, sagging skin and large pores and other problems.

# 8 standard functions



## Small bubble

Breaking through the technical level of traditional small bubble instruments, opening pores, softening keratin, cleansing the skin deeper, improving dullness, and blackheads leaving the skin clean, translucent, and hydrated.



## Bipolar warm wave

Used for facial thermal therapy to allow skin tissue to absorb a large amount of heat energy, combined with new collagen, to smooth and tighten the skin.



## Five polar EMS

Using 20Hz–1000Hz frequency conversion EMS, the face is lifted and tightened, improving fine lines in the eyes, and restoring Eyes expression!



## Bipolar warm wave

Using 20Hz–1000Hz frequency conversion EMS, the face is lifted and tightened, improving fine lines in the eyes, and restoring facial expression



## High Pressure Oxygen Injection Gun

High pressure oxygen injection penetrates deep into the skin layer, making it easier to absorb by narrowing pores.



## Point vibration shovel

The shovel handle utilizes high-frequency vibration mechanical waves to soften the skin's keratinocytes and deeply clean skin dirt.



## Flat head Ultrasound

Using sound wave resonance frequency of 1–3 million vibrations to open pores and allow products to penetrate deep cells.



## Cold import

Quickly conduct skin qualitative analysis and achieve the effect of removing redness, anti allergy, and firming and shaping.

# Specifications —

Power:  $\leq 100W$

Voltage: 100~240V 50Hz /60Hz

Ultrasound: 1.1MHz

Cold: Cold 0°C–6°C

Ultrasound: 28~30KHz

HR Frequency: 0.2–0.5Mhz

Clean Water Mill: Positive pressure 4–5Kg

Negative Pressure 90Kpa , Flow  $\leq 30L$

Display: 12.1 inches, 1024x600

Small bubble product bottles:



Purifying Serum



Hydrating Serum



Revitalizing Serum




Cleansing Liquid



**Micro Diamond Probe: 3 soft heads, 8 hard heads**



# The Effect of Continuous Use



# 9 optional configurations




## Electric Microcrystals

Using 20Hz–1000Hz frequency conversion EMS, four roller massage, facial lifting and firming, improving fine lines in the eyes, and restoring facial radiance!



## Radiant oxygen bubbles

Create an oxygen supply environment, activate cells to enhance absorption, and inject products into deep layers of the skin; Collaborate with fruit acid skincare products.



## Electric microcrystals

The front-end microcrystals or probes impact the surface of the skin back and forth, Generate channels for product import, Effectively promoting product absorption



## LED Colorful Light Mask

Red light promotes collagen regeneration and anti-aging. Blue light sterilization, repairing damaged skin. Green light improves redness and pigment redness. White light accelerates tissue metabolism and improves fine lines.



## Point Head Ultrasound

The front-end microcrystals or probes impact the surface of the skin back and forth. Generate channels for product import and effectively promoting product absorption.



## Diamond Dermabrasion

The principle of diamond particle friction and vacuum suction to remove the aging keratin from the outer layer of the skin epidermis. And promote skin renewal, non-invasive, non irritating, safe and convenient.



## Ozone Ion

Anti inflammatory, antibacterial, soothing and relieving itching, targeting allergic skin, improving skin condition.



## Beauty Spotted Pen

Remove and sweep spots, improve skin radiance.



## High frequency wave

Anti inflammatory, bactericidal, soothing, and relieving itching For allergic skin, Improve skin



**Head Office:** Office 509, Churchill Executive Tower, Al Amal St., Business Bay, Dubai,UAE  
**Phone:** +971 4 569 3474 **Mobile:** +971 50 912 4577 **Email:** [admin@skyline-med.com](mailto:admin@skyline-med.com)